



Disk Filter

Major Applications

Small quantity filtration of various liquids

Features

- Selectable from various filter product types
- Disk filters are packed in a special case for each part number.

Advantages

- Ideal for small volume filtration tests for filter selection
- Easy to control the number of filter media without mishandling

Ordering Information

Filter Cartridge Product Type

S H P

*Please refer to the table on the back.

Micron Rating

- 0 3 0 -

* Please refer to each filter catalog.

Support media Specifications

M

M = No support media
P = With down-stream side support media

Disk Size

4 7

25 = ϕ 25mm
47 = ϕ 47mm

Product Type

DKP

Packaging Code

X

X = 15 sheets / box (*1)
D = 15 sheets / box
E = 20 sheets / box
F = 25 sheets / box
G = 50 sheets / box

*1 Packaging CODE X is
15 sheets/ box =
5 sheets \times 3 cases / box

*2 : There may be a sponge inside the case to prevent the filter media from turning over.
*3 : The sponge will turn yellow when exposed to light, but this does not affect the quality.
*4 : Order is "box" unit.

Filter product type available table

Filter Cartridge Product Type	Support media Specifications	Packaging code for sale	sheets / box
ACES	P	G	50
ACTX	P	G	50
AMPA	M	F	25
AMPX	M	E	20
ASCP *1	M	X	15
ASP *1	M	X	15
CEM	P	F	25
CER	M	G	50
CES	P	G	50
CLP	M	E	20
CTA	M	G	50
CTB	M	G	50
CTY	M	G	50
CTZ	M	G	50
CVP	M	F	25
MAP	M	E	20
MBP	M	F	25
MFB	M	F	25
MGB	M	E	20
MGD	M	F	25
MNN	M	F	25
MPN	M	F	25
MPS	M	F	25
MPT	M	F	25
MPX	M	E	20
NEN	M	F	25
NET	M	F	25
PED	M	D	15
PSG	P	G	50
PSGD *3	P	F	25
SBP *1	M	X	15
SCP *1	M	X	15
SGP *1	M	X	15
SHD *1	M	X	15
SHP *1	M	X	15
SHPH *1	M	X	15
SLP *1	M	X	15
SLPC *1	M	X	15
SLPG *1	M	X	15
SNP *1	M	X	15
STP *1	M	X	15
ZCES	P	G	50
ZCESD *3	P	F	25

*1: We recommend the use of our special disc holders for the subject disc filters.
For details, please contact our sales staff.

*2: Please contact us for the availability of products not listed in the table.

*3: Please use ZCES and PSG for bubble point measurement and bacteria challenge test.

*The contents of the catalog are subject to change without notice.

*The performance data listed in the catalog are Typical values obtained under specific conditions based on our tests.

ROKI GROUP International Pte. Ltd.

6-20-12, Minami-Oi, Shinagawa-ku Tokyo, 140-0013 Japan
TEL: +81-3-5764-1131 FAX: +81-3-5764-0681

www.rokiglobal.com

Manufacturing is based on our Quality Management Systems that meet ISO9001 standards.

Scope

Design, Development, manufacture, and sales of filter cartridges, housings and filtration equipment.



JQA-QMA16323

MS CM009

4th Issue
Disk_210401E